


# PreadyTake, an *in vitro* ready-to-use cell-based model to evaluate potential drug-transporter interactions

Jonatan Cucala, Marta Ollé, Àlex C Estivill, Sheila Guisado and Lourdes Gombau  
ReadyCell SL, Barcelona, Spain (contact: [jucucala@readycell.com](mailto:jucucala@readycell.com))

## INTRODUCTION

Drug transporter proteins may compromise drug permeability across body barriers. Among the large number of transporters, **special attention** has been given to the **Solute Carrier Transporter (SLC)** family because of their role in the **renal and hepatic elimination of drugs**. Understanding the role of these transporters is normally performed *in vitro*, although tools are limited and may not reflect the true impact of a transporter on drug disposition. **PreadyTake** is a family of **ready-to-use HEK293 cell-based models individually expressing hepatic (OATP1B3), and renal (MATE1, OCT2) transporters** that emerges as useful tool to **evaluate potential drug-transporter interactions**.




Localization of the membrane-associated OCT2, MATE1 and OATP1B3 transporters in renal and hepatic cells.

## OBJECTIVES


The main objective of this study is to probe PreadyTake as an *in vitro* ready-to-use model to assess OATP1B3- MATE1- and OCT2- substrates, inhibitors and drug-transporter interactions.

## MATERIALS AND METHODS

### Liquefaction and Shipping Medium Exchange



### Uptake Assay



## RESULTS

Data were in compliance with FDA guidelines on drug-transporter interactions. Reference compounds uptake were at least two-fold that of cell expressing the empty vector (HEK293-MOCK). Furthermore, **absorption decreased by more than 50%** when reference compounds were incubated in the presence of transporter inhibitors.




Figure 1. OATP1B3-mediated valsartan internalization in the absence (Panel A) or presence (Panel B) of cyclosporin A, an OATP1B3 inhibitor.




Figure 2. MATE1-mediated metformin internalization in the absence (Panel A) or presence (Panel B) of quinidine, a MATE1 inhibitor.




Figure 3. OCT2-mediated MPP<sup>+</sup> internalization in the absence (Panel A) or presence (Panel B) of doxepin, an OCT2 inhibitor.

Data indicate that PreadyTake are compliant and useful *in vitro* tools to screen OATP1B3, OCT2 and MATE1-mediated drug-transporter interactions and/or induced drug hepatic/renal toxicity at the early stages of drug development.



[readycell.com](http://readycell.com)  
[reagents@readycell.com](mailto:reagents@readycell.com)  
ReadyCell

**ReadyCell**  
EASING YOUR SCREENING

### REFERENCES:

- Transporter-Mediated Drug Interactions Guidance for Industry, 2020. Food and Drug Administration (FDA). *In Vitro Drug Interaction Studies – Cytochrome P450 Enzyme- and Interactions Guidance for Industry, U.S.*
- Nies AT et al. *Arch Toxicol* 90:1555–1584, 2016.
- Tweedie et al. *Clinical Pharmacology & Therapeutics*, 94(1), 113-125, 2013.
- Hermann Koepsell et al. *Pharm Research*, 24:1227-1251, 2007.